

# CROSSVISION FORTE

"CROSSVISION FORTE" is a sterile solution for ophthalmic use with a physiological pH, based upon a patented technology of cross-linking Hyaluronic Acid (HA) with Urea, therefore enhancing and improving the specific mechanical properties of native HA, resisting to hyaluronidase degradation much longer than native HA, thus with higher and prolonged lubricant and hydrating action than native HA.

The hydrating web structure of cross-linked HA is capable to deliver water for a prolonged period improving the lubricating properties of native HA while boosting its moisturizing effects at the same time it is therefore indicated as an adjuvant in dry-eye syndrome. The cross-linked HA provides to the product lubricating, moisturizing and repairing properties, capable of giving relief to ocular dryness due to dry eye syndrome, irritated or tired eye in case of inflammatory conditions or environmental factors.

The formula was developed to ensure maximum comfort already in the first moments after the application. It is suitable for everyday use, daily and night-time.

It does not contain preservatives and/or molecules with anti-inflammatory and vaso-constrictive actions.

"CROSSVISION FORTE" is compatible with contact lenses.

## INTENDED USE

"CROSSVISION FORTE" are indicated:

- for reinstatement and improvement of tear film stability in dry eye syndrome;
- after refractive surgery;
- to relieve the eyes in cases of reddening, fatigue and discomfort of the eye or the ocular mucosa due to:
  - environmental factors (air conditioned, excessive exposure to solar radiation, wind, smoke, pollution, sea water or pool water and dried salt);
  - dryness of the eyes due to excessive fatigue caused by working on the computer for lengthy periods, study;
  - lacrimal abnormalities secondary to inflammations, chronic blepharitis;
  - lacrimal abnormalities secondary to eye surgery of various nature or eyelid aesthetic surgery;
  - wearing soft or semirigid contact lenses.

## COMPOSITION

Sodium Hyaluronate Crosspolymer Aqueous Solution (of which Sodium Hyaluronate Crosspolymer 0,4%), Trometamol, Hydrochloridric acid 2N, Sodium chloride, Water for Injections. Preservative free.

## INSTRUCTION FOR USE

- Wash and dry your hands thoroughly before use;
- Open the protective aluminum pouch.
  - Detach the single-dose container from the strip;
  - Twist and remove the cap to break the seal;
  - Place the cap taking care to put it on a clean surface;
  - Instill 2 drops in each eye by gently squeezing the single dose container, unless otherwise prescribed;
  - After use immediately close the single-dose container: the content must be used within 12 hours after the first opening.

## WARNINGS AND PRECAUTIONS

- The product is for single patient use;
- The product is for ophthalmic use only;
- Do not use the product if the aluminum pouch is damaged or if the single dose is not perfectly closed;
- Do not use if are visible particles and/or changes in the color of the solution;
- During the application, do not touch your eyes and any other surface with the dropper of the container;
- Do not use after the expiry date indicated on the package and on the single dose;
- In case side effects during use, interrupt the treatment and consult a doctor;
- Do not use in case of known incompatibility towards any of the components of the solution;
- No interactions with drugs or other eye drops are reported; however, leave at least 15 minutes between administrations of any other eye drops;
- Dispose of the product according to your local waste management regulations;
- Keep out of the reach and sight of children.
- Consult your doctor before you use this product if you are pregnant or breast-feeding;
- For external use, do not swallow;
- It is necessary to report any serious incident occurred in relation to the device to the manufacturer and the Competent authority of the Member state in which you are established.

## SIDE EFFECTS

Immediately after instillation, the product may cause a slight temporary blurring of vision due to the viscosity of the solution. It is therefore advisable to wait for a few seconds before any activity that requires clear vision.

## STORAGE CONDITIONS

Store between 2°C and 30°C in a dry place away from heat and direct exposure to sunlight.

The content must be used within 12 hours after the first opening.  
Store the other single dose vials into the original aluminum pouch until needed.

## PACKAGING

From 5 to 30 re-closable single-dose of 0.5 ml.

## VALIDITY

With un-open pack, not exceed the expiry date stated on the box and on the single dose



DT/01/DT/010



**C.O.C. Farmaceutici S.r.l.**

Via Modena, 15 - 40019

Sant'Agata Bolognese (BO) Italia

Last Revision: 02.01 del 10/2022



**C.O.C. Farmaceutici srl**

Via Modena, 15  
40019 - S. Agata Bolognese (BO) ITALIA  
Via Chiesa Sud, 156 C/D/E/F/G  
41016 - Rovereto s/S - Novi di Modena (MO) ITALIA

Cap. Soc. € 1.200.000 I.V.  
C.F. e Iscr. Reg.  
Impr. di Bologna 03541870378  
P.IVA IT 01603491208  
R.E.A. Bologna N. 299432



**Tubilux Pharma S.p.A.**

Via Costarica, 20/22  
00071 Pomezia (Roma)  
ITALIA

cap. soc. € 472.500,00 I.V.  
C.F. / P.IVA 05406661008  
R.E.A. Roma N. 883476