

# SWEETCLEAN SPRAY BABY

## INTRODUCTION

**Sweetclean spray baby** is an isotonic nasal spray without preservatives containing sea minerals from the Gulf Stream, chamomile distilled water and aloe vera gel.

## DESCRIPTION

**Sweetclean spray baby** is a sterile, isotonic, physiological pH saline solution without preservatives, designed to hydrate and protect the nasal mucosa, encouraging fluidification and elimination of nasal secretions.

Sea minerals from the Gulf Stream water help clean the nasal cavity by removing mucous and reducing irritation of the mucosa.

Aloe vera gel, known for its protective and hydrating action, creates a barrier film that protects the nasal mucosa.

Chamomile distilled water works in tandem with the soothing action of Aloe.

**Sweetclean spray baby** contains no preservatives and/or vasoconstrictive substances.

**Sweetclean spray baby** contains no gas propellants.

## TARGET USE

**Sweetclean spray baby** is recommended for the daily nasal cleaning of adults and children (from 6 months) by helping to clear the nose, even when suffering from a cold, and improve breathing. This formulation is useful for hydrating dry nasal mucosa and gently soothing irritation.


## USE INSTRUCTIONS

1. Remove the safety cap (fig. 1);
2. Hold the bottle with three fingers of one hand (fig. 2):
  - hold the spray pump between the index and middle fingers and support the bottom of the bottle with the thumb.
3. When first used, press the dispenser until the spray is regular and even;
4. Insert the tip of the dispenser into the nostril and spray while keeping the head upright and the bottle in vertical position with the thumb on the bottom and the two fingers on the pump (fig. 3);
5. Repeat in the other nostril;
6. Allow the product to work for a few minutes, then, if necessary, blow the nose or wipe away excess mucous;
7. After use, rinse the dispenser with water, dry and replace the safety cap (fig. 4).

Use twice a day, or in accordance with medical advice.

After first use, the bottle contents should be used within 3 months.

Strictly personal use only.



## COMPOSITION

Sea minerals from the Gulf Stream water, Aloe vera gel, Chamomile distilled water, Sodium citrate, Citric acid, Water for injections.

## PACK

The pack contains 1 plastic bottle with sterile contents, a spray dispenser, and no preservatives.

## WARNINGS

Solely for nasal use.

The contents of the spray bottle should be used by one person only.

Do not use the product if the bottle is not perfectly sealed and/or intact.

Do not tamper with or damage the bottle and its mechanisms.

Do not use the product if you are allergic to any of the components in the formulation.

If side effects occur during use, stop using the product and consult a doctor.

Seek advice from a doctor if wishing to take inhaled medicines at the same time as the nasal spray to assess the possibility of incompatibility.

Do not use after the expiry date shown on the pack.

Do not swallow.

Dispose of container correctly after use.

Keep out of reach of children.

Keep dry, away from direct heat and sunlight, and at a temperature between 2°C and 30°C.

## VALIDITY

If unopened, do not use the product after the expiry date shown on the pack and on the bottle label.



**C.O.C. Farmaceutici S.r.l.** Via  
Modena, 15 - 40019 Sant'Agata  
Bolognese (BO) Italy



MEDICAL DEVICE

Last revision: 02.01 of 02/2020



**C.O.C. Farmaceutici srl**

Via Modena, 15  
40019 - S. Agata Bolognese (BO) ITALIA  
Via Chiesa Sud, 156 C/D/E/F/G  
41016 - Rovereto s/S - Novi di Modena (MO) ITALIA

Cap. Soc. € 1.200.000 I.V.  
C.F. e Iscr. Reg.  
Impr. di Bologna 03541870378  
P.IVA IT 01603491208  
R.E.A. Bologna N. 299432



**Tubilux Pharma S.p.A.**

Via Costarica, 20/22  
00071 Pomezia (Roma)  
ITALIA

cap. soc. € 472.500,00 I.V.  
C.F. / P.IVA 05406661008  
R.E.A. Roma N. 883476